



P710LERM36SF



P7102PM36SEF



P710LNQM36SEF



PRODUCT

## Lab Tables

FEATURES

- Available in a variety of sizes and different top options.
- Colored 3mm PVC edge banding options available on HPL and ChemGuard top options.
- 2-1/4" square legs, with 5/16" hanger bolts installed, attached with supplied hardware.
- Aprons are 13/16"D x 4-1/2"H with heavy-duty 13-gauge steel corner braces.
- Solid oak or maple hardwood sealed with an earth-friendly UV finish.
- Equipped with non-skid adjustable glides for precise leveling and rubber leg boots for added protection.
- Choose from plain, drawer or book compartment aprons.
- Book compartment apron tables have either two or four book compartments (depending on size) in the apron measuring 14"W x 13-5/8"D x 3-7/8"H.
- Half and full shelf options: 3/4" thick plywood, veneer finish with edge-banding on all 4 sides.
- Shipped KD with tops attached.
- Weight capacity is 500 lbs evenly distributed

7000 SERIES TABLE CONFIGURATOR

APRON	LEG TYPE	SIZE	TOP TYPE	HPL COLOR (L)	3MM EDGE-BANDING (2/L)	WOOD	HEIGHT	OPTIONS (REQ)
Plain Apron (P)	Wood (7)	48" x 24" (10)	Black T-Mold (1)	Black (B)	Black (B)	Oak (K)	30	Glides (N)
Drawers (D)		48" x 42" (11)	ChemGuard (2)	Charcoal (C)	Charcoal (C)	Maple (M)	34	Leg Braces (L)
Compartment (C)		54" x 30" (13)	Phenolic (4)	Almond (A)	Almond (A)		36	Electrical (E)
		60" x 30" (14)	Maple BB (5)	White (W)	White (W)			Leg Braces, Stretcher (S)
		72" x 30" (15)	Epoxy (6)	Natural Maple (N)	Natural Maple (N)			Leg Braces, Elec (LE)
		48" x 21" (16)	ShopTop (7)	Erasable (E)	Carrot (T)			Leg Braces, Stretcher, Elec (SE)
		48" x 36" (17)	HPL (L)		Lime (L)			Leg Braces, Stretcher, Casters (SC)
		54" x 21" (18)			Aqua (Q)			Leg Braces, Stretcher, Full Shelf (SF)
		54" x 36" (19)			Red (R)			Leg Braces, Stretcher, Half Shelf (SH)
		54" x 24" (20)			Peacock (P)			Leg Braces, Stretcher, Elec, Full Shelf (SEF)
		60" x 21" (21)			Royal (Y)			Leg Braces, Stretcher, Elec, Half Shelf (SEH)
		60" x 36" (22)						Leg Braces, Stretchers, Casters, Full Shelf (SCF)
		72" x 24" (23)						Leg Braces, Stretchers, Casters, Half Shelf (SCH)
		72" x 36" (24)						
		72" x 24" (30)						
	72" x 42" (40)							
	48" Rd (48)							
	60" x 24" (60)							
	54" x 42" (80)							
	60" x 42" (90)							